

PART NUMBERS AND IMAGES

The following dispenser guns can be used with the part numbers assigned.

110007-12



CoMax



Kerr (Original Model)



Caulk



Clinician's Choice

110007-02



Centrix



Venus

110007-06



Kuraray



Vident

110007-07



Centrix Nu-Gun



Vivadent

OTHER ADDENT PRODUCTS



◀ RITE-LITE™
L.E.D Shade Matching Light
PIN 110047



◀ TRIMAX™
Composite Instrument
PIN 110001



◀ MICROLUX/DL™ Kit
with 6 bottles of solution
PIN 110035



◀ MICROLUX™
Fiber-optic mirror
attachment
PIN 110045



◀ MICROLUX™
Transilluminator
PIN 110019



◀ CoMax™
Composite Dispenser
PIN 120000



◀ PowerBrite™
Sufficient for 16 teeth
PIN 111011

CALSET™

COMPOSITE WARMER



Dr. Paul Belvedere says:

"The contemporary heavy bodied materials which are more difficult to express at room temperature become very user friendly with the Calset."

Dr. Frederick A. Ruggenberg says:

"The pre-heating of composite, such as with the Calset, provides significant increase in composite cure and lowers paste viscosity allowing for enhanced composite flow over using the material at room temperature."

Dr. Douglas Terry says:

"As the search for an ideal composite continues, an additional variable has been introduced to direct resin technology- elevated temperature. Elevated temperature in direct placement of composite resin with such instruments as the Calset heater enhances the sculptability and the ease of placement of the composite while also providing an increased conversion and decreased light exposure time for the resin."

AdDent, Inc.

43 Miry Brook Road Danbury, CT 06810

(203) 778 0200 (203) 792 2275

www.AdDent.com info@AdDent.com

P/N 660024 7/2010

